


## Protocol for shRNA construction-I: PCR method

### Preparation of cloning vector:

1. Incubate 3  $\mu\text{g}$  of shRNA cloning vector with 5 units (NEB) of *EcoRI* and 10 units of *AgeI*, (double digestion) in a reaction volume of 100  $\mu\text{l}$  at 37<sup>o</sup>C for overnight (using NEB #4 buffer).
2. Take 5  $\mu\text{l}$  of reaction mixture to check the digested DNA in a regular agarose gel electrophoresis. An example of such restriction analysis is shown in below:



**Figure legend:** 400ng of various shRNA cloning vectors were digested with 2.5 units of *EcoRI*, 5 units of *AgeI* in a reaction volume of 30  $\mu\text{l}$  at 37<sup>o</sup>C for 6 hrs (using NEB #4 buffer), and then 15  $\mu\text{l}$  of mixtures mixed individually with loading buffer were subjected to electrophoresis (0.5% agarose gel containing 0.1 $\mu\text{g/ml}$  of EtBr, run at 50 volts in 0.5X TBE running buffer containing 0.1 $\mu\text{g/ml}$  of EtBr) for 2 hours. As shown in this figure, the 1.9 kb stuffer sequence was completely digested out under this reaction condition, an indication of complete digestion of those vectors.

**Please note that *EcoRI* is easy to trigger star activity, an altered or relaxed specificity of the enzyme (please refer to NEB catalog and technical reference for how to avoid star activity of restriction enzymes).**

3. Resolve restriction enzymes digested shRNA cloning vector by 0.7% agarose gel (mix of 0.35% regular agarose plus 0.35% low melting agarose recommended).
4. Purify DNA from gel according to the instructions of Roche or Qiagen DNA elution kit and elute it with 50  $\mu\text{l}$  of autoclaved 0.1X TE buffer.

## Preparation of insert (PCR product):

1. Design and order a long oligonucleotide as a PCR template for amplifying shRNA sequence you desired to clone (see appendix 1 for how to design shRNA oligonucleotide sequence).
2. Set up PCR reaction:
 

shRNA oligonucleotide (50 $\mu$ M)	1 $\mu$ l
Forward primer (100 $\mu$ M)	2 $\mu$ l
Reverse primer (100 $\mu$ M)	2 $\mu$ l
KAPAHiFi DNA Polymerase	1 $\mu$ l
5X GC buffer	10 $\mu$ l
2mM dNTPs	5 $\mu$ l
d.water final volume to	50 $\mu$ l
3. PCR parameter:
 

Step1 $\rightarrow$ 95 $^{\circ}$ C, 1cycle	5-min
Step2 $\rightarrow$ 95 $^{\circ}$ C	15 Sec
Step3 $\rightarrow$ 65 $^{\circ}$ C	30 Sec
Step4 $\rightarrow$ 72 $^{\circ}$ C	20 Sec
(Repeat step 2 to step 4 for another 2 cycles)	
Step5 $\rightarrow$ 72 $^{\circ}$ C, 1 cycle	5-min
Hold at 20 $^{\circ}$ C	
4. Purify PCR product using MinElute Gel Extraction Kit (Qiagen), and elute it with 80  $\mu$ l of autoclaved 0.1X TE buffer.
5. Digest eluted PCR product with *Bsm*BI and incubate at 55 $^{\circ}$ C in the incubator rather than in water bath for O/N (see appendix 2 for the restriction pattern of *Bsm*BI restriction enzyme).
 

Digestion condition:

Eluted PCR product	80 $\mu$ l
NEB #3 10X buffer	10 $\mu$ l
<i>Bsm</i> BI (10U/ $\mu$ l; NEB)	2 $\mu$ l
d.water final volume to	100 $\mu$ l
6. Purify *Bsm*BI-digested PCR product using MinElute Gel Extraction Kit (Qiagen) and elute it with 20  $\mu$ l of 0.1X TE buffer. (Recovery rate for 55-60 bp DNA is satisfactory if using this kit.).
7. Set up the ligation reaction (use of 2  $\mu$ l vector and 4  $\mu$ l of PCR products, respectively).
8. Transform in *E. coli* ([Stbl III \[Invitrogen\] recommended; Lentivector is hard to transform into DH5 \$\alpha\$](#) ).

**Protocol for shRNA construction-II: annealing method**

1. Design sense (tail with **CCGG** [*AgeI* cohesive end] sequence at the 5' end) and antisense (tail with **AATT** [*EcoRI* cohesive end] sequence at the 5' end) shRNA oligonucleotide as follows:

If consider **cgcatacgacgattctgtgat** as the target sequence, then sense (up-strand sequence of following example) and antisense (low-strand sequence of following example) oligonucleotide, respectively, will be as follows:

<u>target sequence/passenger strand</u>	<u>siRNA sequence/guide strand</u>
<b>ccgg</b> <u>cgcatacgacgattctgtgat</u> <b>ctcgag</b>	<u>atcacagaatcgtcgtatgcggttttt</u>
gcgcatgctgctaagacacta <b>gagctc</b>	<u>tagtgtcttagcagcatacgcaaaaa</u> <b>ttaa</b>

**ctcgag**: loop sequence of shRNA

**Please note that the design is in concert with TRC shRNA library's design.**

2. Order oligonucleotides with 200 nmole scale and OPC or PAGE purification (oligonucleotide produced by Mission Biotech recommended).
3. Dissolve oligonucleotides into 100 µM with autoclaved distilled water.
4. Prepare 10X annealing buffer:
  - 1M K-acetate
  - 0.3M HEPES-KOH pH7.4
  - 20 mM Mg-acetate
5. Set up annealing mixture:
 

Sense oligo	9 µl
Antisense oligo	9 µl
10X annealing buffer	2 µl
6. Anneal mixture by PCR machine using the following parameters:
  - 95°C, 78°C, 74°C, 70°C, 67°C, 63°C, 60°C, 56°C, 53°C,
  - 50°C, 48°C, 46°C, 44°C, 42°C, 40°C, 39°C, 37°C,
  - 36°C, 35°C, 34°C, 33°C, 32°C, 31°C-----5 min in each step
  - 30°C, 28°C, 26°C, 24°C, 22°C, 20°C-----10 min in each step
  - Hold at 4°C
7. Dilute annealing mixture 50x or 100x with ddH<sub>2</sub>O
8. Set up ligation reaction mixture and ligation for O/N:
 

RE-restricted shRNA cloning vector	2 µl
(Prepare vector followed aforementioned protocol)	
Annealed oligonucleotides	2 µl
10X ligation buffer	1 µl
ligase (1 unit/ µl)	1 µl
d.water final volume to	10 µl



**Appendix 5: Terrific Broth 配方** (Please refer to Molecular cloning for detailed)

Tryptone	12 g/L
Yeast extract	24 g/L
Glycerol(100%)	4 ml/L
KH <sub>2</sub> PO <sub>4</sub>	2.31 g/L
K <sub>2</sub> HPO <sub>4</sub>	12.54 g/L

Ampicillin or Carbenicillin final concentration : 100mg / L