


**Protocol of sample preparation for pooled shRNA half-hairpin microarray**



**【PCR amplification】**

**Step 1 – Target/Puro Enrichment**

**(1-1) PCR**

* gDNA (total 130 µg= about ~500x oversampling of 82k library )	X µl
**Puro F primer (100 µM)(5'-TTCACCGAGGGCCTATTTCCCATG)	6.5 µl
**Puro R primer (100 µM) (5'-CGTGAGGAAGAGTTCTTGCAGCTC)	6.5 µl
10X ExTaq Buffuer (TaKaRa)	65 µl
dNTP (2.5 mM each)(TaKaRa)	52 µl
Ex-Taq (Cat.No.KA4101AA TaKaRa)	3.25 µl
DMSO	40 µl
ddW	650-Xµl
<hr/>	
Total	650 µl

> Aliquot 100µl of PCR mixture to 7 wells attached strip caps PCR tube.

**(1-2) PCR profile**

94°C x 5min → [94°C x 30sec → 55°C x 30sec → 72°C x 1.5min ] x 15 cycles → 72°C x 10min

**(1-3) Concentration**

**Qiagen PCR purification Kit (Cat. No. 28106)**

- > Pool 650 µl of PCR product together and expand with 5x volume of Buffer PB (3.25ml).
- > Add 600 µl each to three Qiagen MinElute columns (total 3 columns), spin at full speed for 1 min, discard the flow through
- > Repeat once, until the all PCR product input to column.
- > Wash each column with 750 µl of Buffer PE
- > Elute with 10 µl of ddW for each column, and pool 30 µl together.
- > Store the Target/Puro enriched PCR product at -20<sup>0</sup>C.

**Step2 – Half-hairpin Probe Generation**

**(2-1) TRC-PCR-I**

Target/Puro enriched template	3.25 µl
Cy5 (control)/Cy3 (experiment) labeled F primer (100 µM) (5'-AATGGACTATCATATGCTTACCGTAACTTGAA)	6.5 µl
Universal unlabeled R primer (100 µM) (5'-TGTGGATGAATACTGCCATTTGTCTCGAGGTC)	6.5 µl
10X ExTaq Buffuer	32.5 µl
dNTP (2.5 mM each)	26 µl
Ex-Taq	6.5 µl
ddW	243.75 µl
<hr/>	
Total	350 µl

- > Aliquot 50µl of PCR mixture to 7 wells attached strip caps PCR tube.

**(2-2) PCR profile**

95°C x 5min → [94°C x 30sec → 50°C x 30sec → 72°C x 1min ] x 35 cycles → 72°C x 10min

**(2-3) TRC-PCR-II (without DNA template !!)**

Cy5 (control)/Cy3 (experiment) labeled F primer (100 μM)	6.5 μl
Universal unlabeled R primer (100 μM)	6.5 μl
10X ExTaq Buffuer	32.5 μl
dNTP (2.5 mM each)	26 μl
Ex-Taq	6.5 μl
ddW	243.75 μl
<hr/>	
Total	254 μl

> Add 50μl of PCR mixture to previous tube, total 100ul/well

**(2-4) PCR profile**

95°C x 7min → 55°C x 2min → 72°C x 60min

**(2-5) XhoI digestion**

PCR product	600 μl
10x NEB buffer 4	80 μl
10x NEB BSA	80 μl
XhoI (NEB)	25 μl
ddW	15 μl
<hr/>	
Total	800 μl

> Digest at 37°C overnight, covered with aluminum foils

**(2-6) Concentration**

**Qiagen PCR purification Kit (Cat. No. 28106)**

> Add 5x volume of Buffer PB (4ml) to the 800μl digested PCR product , and mix well.

> Add 750μl each to 3 Qiagen MinElute columns, spin at full speed for 1 min, discard the flow

through .

- > Repeat the previous step.
- > Wash each column with 750µl of Buffer PE, spin at full speed for 1 min, and discard the flow Through.
- > Spin at full speed for another 3 min.
- > Transfer the column to a clean 1.5ml microcentrifuge tube.
- > Elute with 10µl of ddW for each column, and pool 30 µl concentrated PCR product together.

### **(2-7) Gel purification**

- > Add 6µl 6x DNA loading dye to the 30 µl concentrated PCR product
- > Gel purification using a 2% DNA agarose gel (in 1x TAE solution) running with 100V,~120min
- > Cut the lower band (~100bp size) with an attempt to minimize the gel size (Cy5 quenches EtBr).
- > Extract DNA out using Qiagen gel purification kit (don't add 2-propanol!) (Cat. No. 28706), elute with 30 µl of ddW.
- > Check concentration with using Nanodrop (with “Microarray”)
  - ✧ Cy5-probe (control) usually has a concentration of 60~90 ng/µl, Cy5 intensity ~0.6-0.7
  - ✧ Cy3-probe (experiment) is 1.5-fold more concentrated, Cy3 intensity ~0.8-1.0

(Repeat probe preparation if concentration is less than 45 ng/µl or signal intensity less than 0.6)

### **【Note】**

#### **The input amount of genomic DNA for stage I-PCR**

~For example~

82,000 shRNA \* 500 replicates = 41,000,000 copies

41,000,000 copies / (333 copies/ng)= 123,123 ng =~130 µg

(at least  $2.1 \times 10^7$  cells for gDNA extraction)

#### **Primer information**

Provided by C5 Core

101/01/05

The primers used in this protocol are all synthesized by IDT (Integrated DNA Technologies, Inc)

Scale: HPLC purification

Store:  $-20^{\circ}\text{C}$  in dark (especially for Cy3/Cy5-labeled primer!)