

Section II: Lentiviral production

Introduction:

This contains the protocol for lentiviral production in 6 cm plate.

Lentiviral production consists of the following steps:

Day 1	Seed 293T cells
Day 2 (pm)	Transfect cells with 3 plasmids
Day 3 (am)	<i>16 hours post-transfection</i> – change media; replace with BSA-containing media
Day 4 (am)	<i>24 hours after media change</i> – harvest virus (1) ; replace with BSA-containing media
Day 5 (am)	<i>24 hours after harvest (1)</i> – harvest virus (2); combine harvest (1) and (2); discard cells

I. Materials

1. Cell line: **293T** packaging cells <*recommended: healthy; avoid confluent culture, and no continuous subculture more than 10 passages*>
2. Growth media: DMEM+10% heat-inactivated FBS (Hyclone) + 1X Pen/Strep
3. BSA-containing media: DMEM+10% heat-inactivated FBS (Hyclone) + 1% (W/V) BSA < Bovine Serum Albumin (Sigma, A9430)> + 1X Pen/Strep
4. Plates: 6 cm tissue culture plate
5. Transfection-quality plasmid DNA--
 TRC library plasmid: pLKO.1-shRNA vector
 packaging plasmid: pCMV-ΔR8.91 (containing *gag*, *pol* and *rev* genes)
 envelope plasmid: pMD.G (VSV-G expressing plasmid)
6. Transfection reagent: TransIT LT1 (Mirus Bio.) < alternative: FuGene 6 (Roche) >
7. Sterile Polypropylene storage tube

II. Instructions

1. **Day 1:** Seed 293T cells in 6 cm tissue culture plate. (seeding density: 1.6×10^5 cells/ml; seeding volume: 5 mL) and incubate cells in the incubator (37°C , 5% CO₂) until the following afternoon.

Note: [1] The cells are seeded in no Pen/Strep containing-growth media.

[2] Before transferring to a culture incubator, keep cells settled for one hour in hood to reduce uneven distribution.

2. **Day 2:** The cell density should be among 50% to 70% confluence for transfection with TransIT LT1 reagent.

- a. Prepare a mixture of three plasmids:

Reagent	per 6 cm plate
pCMV-ΔR8.91	2.25 µg
pMD.G	0.25 µg
hairpin-PLKO.1	2.5 µg
OPTI-MEM to total volume	250 µl

Note: The cell density at transfection, DNA concentration for transfection and reagent to DNA ratio should be optimized according to the different transfection reagent kit.

b. Dilute transfection reagent:

Reagent	per 6 cm plate
TransIT-LT1	15 µl
OPTI-MEM to total volume	250 µl

Add the TransIT-LT1 reagent dropwise into OPTI-MEM, mix by swirling the tip or gently flicking the tube (do not mix by pipetting), and incubate at RT for 5 min.

- c. Add the plasmid mixture dropwise into the diluted transfection reagent and mix by swirling the tip or gently flicking the tube
 - d. Incubate the transfection mix for 20-30 minutes at room temperature.
 - e. Carefully transfer the transfection mix to the cells. 293T cells would be sensitive to perturbation, therefore take care not to dislodge the cells from the plate.
 - f. Incubate cells in the incubator (37 °C, 5% CO₂) for 16 hours.
3. **Day 3:** Change media to remove the transfection reagent and replace with **5 mL** BSA-containing media per plate. Incubate cells for 24 hr (37 °C, 5% CO₂).
4. **Day 4:** Harvest media containing lentivirus at 40 hours post-transfection. Transfer media to a polypropylene storage tube and store at 4 °C. Replace with **5 mL** BSA- containing media for further harvest.
5. **Day 5:** Harvest media containing lentivirus at 64 hr post-transfection and then discard the cells. Pool the viral harvests and spin the media at 1250 rpm for 5 minutes to pellet the cells that were collected during harvest. Transfer the supernatant to a polypropylene storage tube and then aliquot to smaller storage tubes to reduce the numbers of freeze/thaw cycles.
Note: Virus may be stored at 4 °C for a few days but should be stored at -80 °C for long-term storage.

Version Notes:

BSA-containing media for viral harvests: By comparison with growth media, viral harvest media containing 1 g/100 mL supplemental BSA can improve virus yield by ~2 fold.