

Section III: Lentiviral Infection**Introduction:**

This section contains protocols for large scale (6 cm plates) and high throughput (96-well and 384-well) lentiviral infections to achieve stable shRNA-mediated target gene knockdown.

Lentiviral infection consists of the following steps:

Day 1	Seed cells
Day 2	Add lentivirus to cells in growth media containing polybrene (optional for 96-well infections: centrifuge cells to promote infection)
Day 2-3	Remove media and replace with fresh growth media
Day 3+	(optional) Select for infected cells with media containing puromycin
Day 5+	Assay infected cells

All lentiviral procedures should be carried out in accordance with biosafety requirements of the host institution.

Part 1: Lentiviral infection in 6 cm plates**I. Materials**

6 cm tissue culture plates (appropriate for cell-based assay)

Human or mouse cell line and appropriate growth media

Reagents required for cell-based assay

Polybrene (Hexadimethrine bromide; Sigma #H9268) or Protamine sulfate (MP Biomedicals #194729)

(Optional) Puromycin Dihydrochloride (Sigma #P8833)

II. Instructions**A. Optimization of lentiviral infection**

Lentiviral infections should be optimized for each cell line and cell-based assay. For example, the following parameters should be tested before starting large-scale infections to determine the optimal conditions for a given experiment:

Cell seeding density

Amount of lentivirus

Puromycin concentration

Timecourse

B. Infection protocol

1. Seed cells at appropriate density in 5 mL media in 6 cm plates.
 - a. Adherent cells: seed 1 day prior to infection. Incubate overnight (37 °C, 5% CO₂).
 - b. Suspension cells: seed at the day of infection in media containing polybrene* (see table in step 2a).

2. Add virus to cells:

- a. (Adherent cells) Remove growth media and add fresh media containing polybrene* (see table). Alternatively, remove a portion of the growth media and supplement with media containing polybrene. Adjust volumes and polybrene concentration to achieve the correct final polybrene concentration.

Reagent	Per 6 cm plate
Media containing polybrene*	to 5 mL
Final polybrene concentration	8 $\mu\text{g/mL}$

Note: Protamine sulfate may be substituted if polybrene is toxic to cells.

b. Add virus to cells.

3. Incubate cells overnight (37 °C, 5% CO₂).

Note: If polybrene or protamine sulfate brings toxicity to cells, then remove media and replace with fresh growth media at infection day.

4. Change media at 24 hours post-infection. Remove media and replace with 5 mL fresh growth media. If puromycin selection is desired, use fresh growth media containing puromycin.

Note: Puromycin concentration should be optimized for each cell line; typical concentrations range from 2-5 $\mu\text{g/mL}$.

5. Incubate cells (37 °C, 5% CO₂), replacing growth media (with puromycin, if desired) as needed every few days. Incubation periods are highly dependent on the post-infection assay. Puromycin selection requires at least 48 hours. The following recommendations are general guidelines only, and should be optimized for a given cell line and assay.

6. Assay infected cells.

Post-infection assay	Incubation time post-infection	Incubation time with puromycin selection
mRNA knockdown (qPCR)	3+ days	2+ days
Protein knockdown (Western)	4+ days	3+ days
Phenotypic assay	4+ days	3+ days

Part 2: Lentiviral infection in 96-well or 384-well plates (high throughput)

I. Materials

96-well or 384-well tissue culture plates

Human or mouse cell line and appropriate growth media

Polybrene (Hexadimethrine bromide; Sigma H 9268) or Protamine sulfate (MP Biomedicals #194729)

(Optional) Puromycin Dihydrochloride (Sigma #P8833)

II. Instructions

A. Optimization of lentiviral infection

Lentiviral infections should be optimized for each cell line and cell-based assay. For example, the following parameters should be tested before starting large-scale infections to determine the optimal conditions for a given experiment:

- Cell seeding density
- Amount of lentivirus
- Puromycin concentration
- Timecourse

B. Infection protocol

1. Seed cells at appropriate density in 96-well (100 μ L per well) or 384-well (50 μ L per well) tissue culture plates.

a. Adherent cells: seed 1 day prior to infection. Allow seeded plates to sit undisturbed at room temperature for 1 hour before transferring to a tissue culture incubator overnight (37 °C, 5% CO₂).

Note: allowing cells to settle at room temperature can reduce uneven distribution of cells.

b. Suspension cells: seed day of infection in media containing polybrene* (see table in step 2a).

2. Add virus to cells.

a. (Adherent cells) Remove growth media and add fresh media containing polybrene* (see table). Alternatively, remove a portion of the growth media and supplement with media containing polybrene to achieve a final polybrene concentration of 8 μ g/mL (following addition of virus).

Reagent	per well, 96-well plate	per well, 384-well plate
Media containing polybrene*	to 100 μ L	to 50 μ L
Final polybrene concentration	Reagent	8 μ g/mL
Virus (added in step 2b)		
High MOI	5 to 10 μ L	2 to 5 μ L
Low MOI**	\leq 1 to 3 μ L	N.D.

* *Protamine sulfate may be substituted if polybrene is toxic to cells.*

** *Low MOI infections may require dilution of virus stock prior to addition to cells.*

b. Add virus to cells (see table in step 2a).

Note: The indicated range of viral volume for high and low MOI infections assume typical viral yields from the 96-well viral preparation method described in Section II.

3. Option 1: Spin infection

- a. Spin cells at 2250 rpm (~1100Xg) in plate for 30 minutes at 37 °C. Centrifugation can improve viral infection and decreases the length of exposure of cells to polybrene and virus.
Note: Centrifugation is not recommended for 6-well plates or larger, as cells may not be fully covered with media during the spin.
- b. (Optional) Change media immediately following spin infection. Remove media and replace with 100 µL (96-well plates) or 50 µL (384 well plates) fresh growth media.
- c. Incubate cells overnight (37 °C, 5% CO₂).
- d. If puromycin selection is desired, remove media 24 hours post-infection and replace with 100 µL (96-well plates) or 50 µL (384 well plates) fresh growth media containing puromycin.
Note: Puromycin concentration should be optimized for each cell line; typical concentrations range from 2-5 µg/mL.

Option 2: No-spin infection

- a. Incubate cells overnight (37 °C, 5% CO₂).
 - b. Change media 24 hours post-infection. Remove media and replace with 100 µL (96-well plates) or 50 µL (384 well plates) fresh growth media. If puromycin selection is desired, use fresh growth media containing puromycin.
Note: Puromycin concentration should be optimized for each cell line; typical concentrations range from 2-5 µg/mL.
4. Incubate cells, replacing growth media (with puromycin, if desired) as needed every few days. Incubation periods are highly dependent on the post-infection assay. Puromycin selection requires at least 48 hours. The following recommendations are general guidelines only, and should be optimized for a given cell line and assay.

Post-infection assay	Incubation time post-infection	Incubation time with puromycin selection
Viral titer (Puromycin selection/cell viability)	3+ days	2+ days
mRNA knockdown (qPCR)	3+ days	2+ days
Protein knockdown (Western)	4+ days	3+ days
Phenotypic assay	4+ days	3+ days

5. Assay infected cells.